

A Quadro Engineering White Paper



Quadro Engineering Corp.  
613 Colby Drive  
Waterloo, Ontario  
N2V 1A1  
Canada  
519-884-9660  
quadrosales@idexcorp.com  
[www.quadro.com](http://www.quadro.com)

## Improving API Particle Size Distribution Targets Is There a Better Way?

*By Wilf Sangüesa, P. Eng.*  
Product Manager - Pharma Segment

## Improving API Particle Size Distribution Targets - Is There a Better Way?

Milling is an integral step in API production; however, not all size reduction methodologies offer the same level of narrow particle size distribution spans, nor an easy, repeatable and accurate ability to control particle size range outputs. As companies strive for better manufacturing efficiencies – selection of the most suitable equipment for the application, needs to be paired with options offering the greatest return-on-investment. This Paper explores the background of API milling from a technology perspective – reviewing traditional comminuting methodologies, and exploring new, innovative breakthroughs in conical milling size reduction.

### Contents

API Milling Trends.....	1
Size Reduction Equipment Alternatives.....	2
Containment.....	3
Conclusion.....	4
References.....	5

### **Active Pharmaceutical Ingredients (API) Milling Trends:**

It is well known, and widely documented, that API production is one of the fastest growing segments within the pharmaceutical industry; with manufacturing production expected to grow at a 7.9% CAGR through 2016.<sup>[1]</sup>

Partly fuelled by the world’s increased need for “secondary care” small molecules – integral in prescription drugs issued by specialists such as oncologists<sup>[2]</sup>, and partly driven by a concurrent increased demand of well-established drugs in both mature and emerging markets<sup>[3]</sup>; API manufacturers have been challenged to look for the most efficient processing equipment available, without sacrificing product quality or operator safety.

The most appropriate equipment selection begins in the early stages of a Drug Discovery Process – usually between pre-clinical and Phase II<sup>[4]</sup>, before submitting a New Drug Application (NDA). It is at these stages that formulators consider a vast number of ingredient properties and characteristics in their decision as to how best to commercially produce and launch their product.

After what is often a costly and lengthy new molecule development period, an optimized API particle size distribution target is reached - where the new drug will offer the best bioavailability, stability, granule properties and technical feasibility.<sup>[5] [6]</sup> It is at this stage, that powder milling or comminuting technologies are investigated and selected according to their ability to duplicate lab results on a larger scale.<sup>[7]</sup>

Over the last century, powder milling equipment has progressed to meet the changing “needs” of APIs, growing safety requirements (of both product and operators), and increasing production volumes - to name but a few of the many drivers that the API industry has seen. However, accelerated changes in the last few decades are trending towards smaller and smaller particle size distributions as potencies increase, and drug delivery methods change.<sup>[8]</sup>

Most milling technologies use varying forms of impact, grinding, crushing or attrition to convert energy into size reduction – whether generated by gyrating impellers or rotors, or by pneumatic means. Since the amount of size reduction is directly proportional to the “energy” being generated within the milling chamber<sup>[10]</sup>, the right equipment selection becomes increasingly important (Figure 1).

Particle sizing, homogeneous dispersion and calibration are of course key functions of a milling device; however, in addition to size reduction steps, API manufacturing requires attention to equipment design, industrial hygiene and safety.<sup>[9]</sup>

### Size Reduction Equipment Alternatives:

Traditionally, API size reduction technologies are selected based on a number of product-specific, and equipment compliant parameters. One of these parameters typically being attainable PSDs .

For comparison purposes, this Paper groups API comminuting technologies into five PSD targets, ranging from “coarser” to “finer” distributions. Oscillators are sometimes utilized at the coarser end of the distribution spectrum (d50s >300 μm), followed by conical mills (typical d50s in the 150-300 μm range), then by hammermills and pin mills in the mid-to-fine-range (d50s ≈ 75-150 μm and 15-150 μm respectively), and finally jet mills for finer PSDs (common d50s < 15 μm).<sup>[11]</sup>

Focusing on the mid-to-fine-range PSD target (d50s ≈ 15-150 μm), hammermills and pin mills have provided adequate results over the last several decades, albeit at times at a trade-off in excess product losses, overheating and/or plugging of internal components. For example, pin mills can lead to very extensive cleaning cycles due




Figure 2: PSD Comparison Between Traditional Conical Mills and High Energy Conical Mills




Figure 1: Relationship Between Particle Size Distribution vs. Amount of Energy Imparted

to their propensity to pin plugging<sup>[12]</sup>, which in turn results in higher percentages of product waste. Additional challenges of these traditional mills include the amount of milled product outside the desired particle size distribution or the difficulty to easily control and shift PSD targets to accommodate varying API characteristics. Furthermore, when looking at pin mills in particular, the need for air handling systems adds to the complexity of the process; with increased product losses, as well as higher production and maintenance costs.

To address some of these processing challenges, recent innovations in conical milling technology have resulted in the ability to produce finer API particle sizes (capable of shifting PSDs targets by 30-70%), taking cone milling performance from d50s of 150-300 μm, down to a sub 150 μm range (Figure 2).

These new, high energy conical mills, not only lower the PSD threshold of conventional cone mills, but also address one of Pharma’s primary concerns: the production of narrower (tighter) particle size distribution spans, thus improving on-target milling by 30-60% (Figure 3).

Furthermore, these high energy conical mills, provide the ability to control and shift PSDs, allowing manufacturers to maximize the amount of milled product within the desired target – reducing (and potentially eliminating) the need for screening over/undersized particles; resulting in less waste (Figure 4).

**Containment:**


Although it is not within the scope of this Paper to cover in detail product and operator safety/containment trends, they nevertheless continue to be one of Pharma’s leading priorities.

Most API milling technologies address this industry pull with a variety of engineered solutions such as pneumatic handling equipment or product collectors for less stringent conditions; to cross-flow booths and isolators for more rigorous applications. Unfortunately, many of these solutions require costly capital investments and result in continually high operating costs. Finding ways to reduce the use of ancillary equipment– while still ensuring containment OEL levels below 1  $\mu\text{m}^3$  adds another variable or criteria to be considered during milling equipment selection.


**Conclusion:**

API producers live in a global environment, steered by the strictest of manufacturing guidelines, quality and safety requirements; and under profitability pressure from increased competition.<sup>[13]</sup>

Over the decades, milling technologies have improved and adapted to changing standards and customer needs – with varying levels of strengths and weaknesses. Choosing the right or optimum milling solution for each API size reduction application can be a critical and laborious balancing act.



**Figure 3:** High Energy Cone Mill Capable of Narrower Particle Size Distributions



**Figure 4:** Controllability and PSD Shift Capabilities of High Energy Cone Mills



**Figure 4:** In-line Tri-Clamp Designed High Energy Cone Mill

On one hand, the properties and characteristics of APIs must be retained; on the other hand, production efficiencies must be maximized.

API manufacturers have an increased number of size reduction equipment choices<sup>[14]</sup>, including new developments in conical mills. When taking into account efficiency savings from tighter PSD spans (less out-of-spec waste), the ability to control/shift PSD targets (closer to the objective) while still ensuring sub 1  $\mu\text{m}^3$  OEL levels – finding a milling solution that delivers “everything” may appear daunting at first. However, new size reduction innovations enable impactful productivity improvements, enviable return-on-investment and an overall improvement in API particle size distribution targets.

*Wilf Sangüesa is Product Manager – Pharma Segment of Quadro Engineering Corp. (Waterloo, Ontario, Canada), and The Fitzpatrick Co. (Elmhurst, Illinois, USA). Wilf holds a degree in Mechanical Engineering from the University of Western Ontario (London, Ontario) and is licensed as a Professional Engineer in the Province of Ontario.*

*Quadro Engineering, and The Fitzpatrick Co. are IDEX Material Processing Technologies Units (Elmhurst, IL), and have been trusted partners of the world's top pharmaceutical, chemical and food ingredient processing customers for decades. Specializing in solid dosage and API milling, roller compaction and security screening, Quadro and Fitzpatrick are known globally for their quality services, application support and leading technology brands, including the Quadro® Comil®, Fitzpatrick's Fitzmill™ and Chilsonator® products. The Fitzpatrick Company (Elmhurst, IL) was founded in 1934, and Quadro Engineering Corp. (Waterloo, Ontario) in 1976. Visit [www.quadrocomil.com](http://www.quadrocomil.com), [www.fitzmill.com](http://www.fitzmill.com), or contact us via e-mail [quadrosales@idexcorp.com](mailto:quadrosales@idexcorp.com) or at (519) 884-9660 for further information.*

---

**References:**

- 1) Active Pharmaceutical Ingredient (API) Market Trends, Competitive Landscape and Global Forecasts (2011-2016); Markets.com
- 2) Datamonitor, PharmaVitae Explorer, January 2012; company-reported information
- 3) MarketLine, Global Generics, August 2013; industry profile
- 4) GAO Report to Congress “New Drug Development, November 2006”
- 5) Pharmaceutical Preformulation and Formulation: A Practical Guide from Candidate Drug Selection to Commercial Dosage Form By Mark Gibson; Published by Informa Health Care, 2001
- 6) Bioavailability and Granule Properties: Sunil S. Jambhekar, South University School of Pharmacy, Savannah, GA, USA ; Handbook of Pharmaceutical Technology, Chapter 19; Taylor and Francis, 2005
- 7) Sizing and Granulation; Gurvinder Singh Rekl, Murall V. Vuppala; Handbook of Pharmaceutical Granulation Technology; Chapter 13; published by Marcel Dekker, 1997
- 8) R & D of Milling Technology in Pharmaceutical Industry; FUKUNAKA TADASHI (Banyu Pharm. Co., Ltd.) Journal of the Society of Powder Technology, Japan; Vol. 40; No. 9; Pg 655-663 (2003)
- 9) Milling of Active Pharmaceutical Ingredients; Elizabeth S. Fisher Merck and Co., Inc., Rahway, New Jersey, U.S.A: Encyclopedia of Pharmaceutical Technology ; October 2006
- 10) Impact milling of pharmaceutical agglomerates in the wet and dry states; Schenck, L.R., Plank, R.V., International Journal of Pharmaceutics 348 (1-2): 18–26 (2007)
- 11) Handbook of Pharmaceutical Granulation Technology, Third Edition (Drugs and the Pharmaceutical Sciences), Dilip M. Parikh, 2009
- 12) “Milling Experiments with Prototypes of The Quadro Fine Grind F-10 Mill”. G. Samburski, Y. Ovardia; CHoPS Conference Italy, Aug 2006 Teva pharmaceutical industries, Israel
- 13) Global Generic Pharmaceutical Market - Qualitative and Quantitative Analysis; Aiswariya Chibambam, Frost & Sullivan; Pharma Tech – Ahmedabad 2013
- 14) Pharmaceutical Dosage Forms: Tablets, Volume 1, Third Edition; Larry L. Augsburger, Stephen W. Hoag, Taylor & Francis 2008